

# AUTOGRAF<sup>®</sup>

WARM SHAPE MEMORY STAPLE



*Body heat activated*



# AUTOGRAF®

The warm shape memory AUTOGRAF® staples are designed for osteosynthesis in hand and fore-foot surgery, when a compression is needed. The large range of sizes fits varying anatomy.

## Advantages

- ✓ Efficient and long lasting compression,
- ✓ No heating required: the implant is activated by body temperature.

## Warm Shape Memory




Freezer



Body temperature

## Ancillaries



## References

| References | Width (mm) | Legs Length (mm) | Thickness (mm) | Drill Diameter (mm) |
|------------|------------|------------------|----------------|---------------------|
| ATM10      | 10         | 10-10            | 1,5x1,5        | 2                   |
| ATB10-15   | 10         | 15-13            | 1,5x1,5        | 2                   |
| ATB10-17   | 10         | 17-15            | 1,5x1,5        | 2                   |
| ATB10-19   | 10         | 19-16            | 1,5x1,5        | 2                   |
| ATM12      | 12         | 10-10            | 1,5x1,5        | 2                   |
| ATB12-15   | 12         | 15-13            | 1,5x1,5        | 2                   |
| ATB12-17   | 12         | 17-15            | 1,5x1,5        | 2                   |
| ATM15      | 15         | 10-10            | 1,5x1,5        | 2                   |
| ATM15L     | 15         | 15-15            | 1,5x1,5        | 2                   |
| ATM17      | 17         | 10-10            | 1,5x1,5        | 2                   |



Campus de Ker Lann  
Rue Blaise Pascal  
35170 Bruz - FRANCE  
Tél +33(0)2 99 05 59 69  
Fax +33(0)2 99 05 95 62  
info@memometal.com  
www.memometal.com  
Patented Implant

CE  
0120